

**An anti-backflow kit be used with a chiller system in the following situations:**

- **Installation Below MRI Equipment:** If the chiller is installed below the level of the processing equipment, an anti-backflow kit prevents any backflow through the chiller when it is turned off<sup>1</sup>.
- **Overhead Piping:** When there is overhead piping **≥32 feet** that could cause a backflow through the chiller when the unit is turned off.

**GE Part#:**

- **E8915LA** - KKT Anti-Backflow kit option for cBoxX 40-100 chillers
- **E8915LB** - KKT Anti-Backflow kit option for cBoxX 120 – 200 chillers

**EXAMPLE INLUSTRATION:** CHILLER Installation Below MRI Equipment

